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Total Non-Reactive/Negative

Total Invalid (or Indeterminate)

R = Reactive, NR = Nonreactive, and I = Invalid

P = Positive, N = Negative and I = Indeterminate

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**PAGE TOTAL**

Total Tested  
 Total Reactive/Positive  
 Total Non-Reactive/Negative  
 Total Invalid (or Indeterminate)

	↓	↓	↓	↓	↓
20	20		20		
4	3		3	3	3
16	17		16		1
0			1		

R = Reactive, NR = Non-reactive, and I = Invalid

P = Positive, N = Negative and I = Indeterminate

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	Test 1	Test 2	% Agreement	EQA	Final status
# Negative	945		100% (assumed)		945
# Positive	55	Pos 53 Neg 2	=53/55 96.4%	Pos 53 Neg 2	Pos 53 Neg 2
Total tested	1000	1000	=998/1000 99.8%		1000

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	<b>Proficiency Testing Program using Dried Tube Specimens (DTS)</b>	
	<b>Doc. No.:</b>	<b>Rev. No.: 00</b>
	<b>Effective Date:</b> August 20, 2007	<b>Page:</b> 1 of 4
<b>Written by:</b> Hetal Patel	<b>Revised:</b> Sept 23, 2008 Bharat Parekh	

## 1.0 TITLE

Proficiency Testing Program using Dried Tube Specimens (DTS)

## 2.0 PURPOSE

The purpose of this procedure is to provide guidance to setup proficiency testing (PT) program using Dried Tube Specimens (DTS) to ensure the quality of HIV testing in PEPFAR countries.

## 3.0 EQUIPMENT

NA

## 4.0 SUPPLIES

- 2.0 mL conical bottom sarstedt tubes (catalog # 72-694-007; [www.sarstedt.com](http://www.sarstedt.com))
- Green food coloring dye (Kroger brand)
- Pipettes those are capable of multi-dispensing (Rainin EDP3 pipette, catalog # E3-20).
- Pipette tips
- Disposable transfer pipettes (Samco Transfer pipet catalog number 282)
- Freezer boxes
- Tube racks

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## 6.5 Preparation of PT buffer (PBS/Tween-20)

- 6.5.1 Phosphate buffered saline with tween 20 pouches can be commercially purchased from Sigma (catalog # 3563).
- 6.5.2 Dissolve one foil pouch of phosphate buffered saline with tween 20, pH 7.4 in 1 L of deionized water. It will yield 0.01 M PBS; NaCl 0.138 M; KCl 0.0027 M; Tween 20 – 0.05%; pH 7.4.
- 6.5.3 Filter it through 0.2  $\mu$  filter flask.
- 6.5.4 Prepare 1.8 ml aliquots in pre-labeled 2 ml screw capped tubes.
- 6.5.5 The label on tube should include the following:
  - 6.5.5.1 Identify the tube as “PT buffer”
  - 6.5.5.2 Set an expiration date of 1 year after you prepare

## 6.6 Preparation of DTS

- 6.6.1 Create a panel of 6 samples from the characterized specimens with a combination of negative and positive reactivity for HIV.
- 6.6.2 Carefully blind the panel assigning a new ID to each of the 6 member panel. For example DTS-A1 to DTS-A6. Ensure there is a link of the original ID and new ID.
- 6.6.3 Label each tube with appropriate new ID.
- 6.6.4 Depending on the number of labs enrolled in the PT program, prepare 10 to 20 extra sets.

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

# Pilotaje de DTS en Kenya

## Número de sitios CT = 24

Table 2. Number of testing sites reporting screening results

Specimen Number	Expected Results	Positive	Negative	Invalid	No interpretation
A1	Positive	24	0	0	0
A2	Positive	24	0	0	0
A3	Negative	0	24	0	0
A4	Negative	0	24	0	0
A5	Positive	24	0	0	0
A6	Negative	0	23	1	0

# Pilotaje de DTS en Kenya

## Conclusion, Recommendation and Way forward

- The sites that participated in PT DTS pilot performance were good.
- We strongly feel that PT for HIV testing using DTS has more merits than repeat testing which is currently ongoing.
  - easy to prepare
  - specimen stable at high temperatures
  - easy storage and transport
  - Short turnaround time for feedback
- Based on outcome of the pilot PT, we strongly recommend an elaborate PT roll-out plan to all HIV testing sites
- PT using DTS should supplement repeat testing. Retesting should be scaled down to new sites and poorly performing old sites.

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