

Regulatory Issues for Rapid Tests in the U.S.

2nd International Workshop

Accessible and Quality Assured In Vitro
Diagnostic Tests for Public Health Programs

November 19, 2012

Elliot P. Cowan, Ph.D.


Disclaimer

This presentation reflects the views of the presenter and should not be construed to represent FDA's views or policies.


Goal of this Talk

- To show how public health need can drive the regulation of HIV tests
 - Professional rapid HIV tests
 - Over-the-counter home-use HIV tests

Evolution of HIV Testing



Evolution of HIV Testing



The Gap: 1990s

- 25% of Americans living with HIV did know they are infected
- 37% first discovered they had HIV ≤ 1 year before an AIDS diagnosis
- 31% who tested positive at CDC-funded HIV testing sites do not return for results
 - Testing required two visits, with results in 1-2 weeks


Filling the Gap: Rapid HIV Tests

- Only one visit required to get a test result
- Relatively non-invasive
- Tests can go to where the people are
 - Not lab or clinic-based
 - Opportunity for non-traditional testing sites

Bringing Rapid HIV Tests to Market

- FDA simplified clinical trial requirements
- FDA identified rational performance expectations
 - Sensitivity and specificity $\geq 98\%$ as the lower bound of the 95% confidence interval

Evolution of HIV Testing



The Gap Persists: HIV in the US (CDC)

Number of Americans newly diagnosed with HIV infection each year



~50,000

Number of HIV-infected people in the US



1.2 MILLION

HIV-infected people in the US who do not know they are infected



1 in 5

Filling the Gap: The Sequel

- Why does the gap exist (barriers to testing)?
 - Stigma
 - Inconvenience
 - Denial
 - Etc.
- Would an over-the-counter rapid HIV test help to bridge the gap?
- Why now?
 - HIV infection perceived as treatable

Over-The-Counter (OTC) Home-Use Test Kits

- An over-the-counter test kit is one that is available for purchase by a consumer for the purpose of self-testing
 - No referral from a medical professional (no prescription needed)
 - Many purchase options
 - Internet purchase
 - Pharmacy/Drug Store/Grocery Store
 - Vending machine

Home-Use Test Kits

- Features
 - Specimen is collected by the test subject
 - Test is performed and interpreted by the test subject
 - No trained operator
 - Results obtained without personal counseling
 - Counseling and medical referral available by phone
 - Screening result needs to be confirmed by an additional professional test

Important Concepts

- Over-the counter tests are COMPLETE test systems that are not only the tests themselves but the entire testing process
 - Design test system for self-testing
 - Validate test performance in the hands of the intended users
- Lower sensitivity and specificity expected in the hands of untrained users compared to professional users.

Expected Performance

- Discussion in public meetings
- FDA Blood Products Advisory Committee recommended sensitivity and specificity $\geq 95\%$ as the lower bound of the 95% confidence interval (vs. 98% for professional use rapid HIV tests)

OraQuick® In-Home HIV Test System




Professional Test Performance

| | Performance of the OraQuick ADVANCE [®] Rapid HIV-1/2 Antibody Test (2-sided 95% CI) | BPAC Minimum Recommended Performance |
|-------------|-----------------------------------------------------------------------------------------------|--------------------------------------|
| Specificity | 99.8% (99.6 - 99.9%) | 98% (lower bound of the 95% CI) |
| Sensitivity | 99.3% (98.4 - 99.7%) | 98% (lower bound of the 95% CI) |

Home-Use Test Performance

| | Performance of the OraQuick® In-Home HIV Test Kit (2-sided 95% CI) | BPAC Minimum Recommended Performance |
|-------------|--------------------------------------------------------------------|--------------------------------------|
| Specificity | 99.98% (99.90 - 100%) | 95% (lower bound of the 95% CI) |
| Sensitivity | 92.98% (86.64 - 96.92%) | 95% (lower bound of the 95% CI) |

Home-Use Test Performance as a Measure of Safety and Effectiveness



Risk Analysis

(FDA)

- Estimated test results for numbers of individuals projected to be tested annually who would not otherwise be tested
 - True positive, false negative
 - True negative, false positive
- Estimated net transmissions averted
- Impact of switching from professional testing to self-testing
- Impact of who will use the test
- Public health implications and individual health implications

FDA Projected Outcomes of Testing with the OraQuick[®] In-Home HIV Test in the 1st Year

(Based on point estimates of sensitivity and specificity)

| | |
|-----------------------|---------------|
| True Positive | 45,000 |
| False Negative | 3,800 |

| | |
|-----------------------|------------------|
| True Negative | 2,700,000 |
| False Positive | 1,100 |


13 TP : 1 FN (vs. 128:1)

↑
Professional
use test

↓
3,750 TN :1 FP (vs. 461:1)

Projected Net Transmissions Averted

(Combined Professional and Home Testing)




Source: FDA presentation at May 15, 2012 BPAC meeting

Summary of FDA Assessment

- Based on a risk assessment model, FDA projected a net public health benefit to the OraQuick[®] In-Home HIV Test
 - Net increase in number of HIV infections newly identified in the first year
 - Net transmissions averted
- Individual risk remains in the form of increased numbers of false negative results

Evolution of HIV Testing



Messages that Came with Approval

- A positive result with this test does not mean that you are definitely infected with HIV, but rather that additional testing should be done in a medical setting.
- A negative result with this test does not mean that you are definitely not infected with HIV, particularly when exposure may have been within the previous three months.
- Resting is recommended if you test negative and continue to engage in behavior that puts you at risk for HIV infection.
- A negative result does not imply it is safe to engage in risk behavior for HIV infection.






Some Things to Think About

- What is the need?
- What tools are available to meet that need?
- What are the gaps left by existing tools?
- How can those gaps be filled?
- How is the system monitored to determine if the right choices were made?
- Do the benefits outweigh the risks?



The Power of Sharing (or Emoticommon Sense)



Let the Public Health Need Drive the Process

